

## Pipetting Secret Code - KEY

	1	2	3	4	5	6	7	8	9	10	11	12
A												O
B											O	
C		O	O		O	O	O	O		O	O	
D	O			O			O		O		O	
E	O			O		O				O	O	
F	O	O	O	O	O	O	O				O	
G												O
H												

Pipet the following amounts in the designated assay plate cells:										
Cell	Yellow	Red	Blue	Green		Cell	Yellow	Red	Blue	Green
D4	0	0	350 µL	0		C7	150 µL	25 µL	25 µL	150 µL
E4	0	0	350 µL	0		D7	150 µL	25 µL	25 µL	150 µL
F4	0	0	350 µL	0		F7	150 µL	25 µL	25 µL	150 µL
C8	0.2 mL	0	50 µL	100 µL		C6	0.1 mL	50 µL	50 µL	150 µL
F8	0.2 mL	0	50 µL	100 µL		E6	0.1 mL	50 µL	50 µL	150 µL
A12	325 µL	25 µL	0	0		F6	0.1 mL	50 µL	50 µL	150 µL
G12	325 µL	25 µL	0	0		C5	50 µL	50 µL	50 µL	200 µL
B11	225 µL	100 µL	0	25 µL		F5	50 µL	50 µL	50 µL	200 µL
C11	225 µL	100 µL	0	25 µL		C3	0	25 µL	300 µL	25 µL
D11	225 µL	100 µL	0	25 µL		F3	0	25 µL	300 µL	25 µL
E11	225 µL	100 µL	0	25 µL		C2	0	75 µL	250 µL	25 µL
F11	225 µL	100 µL	0	25 µL		F2	0	75 µL	250 µL	25 µL
C10	75 µL	0.25 mL	0	25 µL		D1	0	100 µL	250 µL	0
E10	75 µL	0.25 mL	0	25 µL		E1	0	100 µL	250 µL	0
D9	0	350 µL	0	0		F1	0	100 µL	250 µL	0